

60*90mm



QUICK GUIDE

C2-Instructions for Use

Smart Cooling Companion

Insulin Pen Cooler Travel Bottle



FREEZE COOLING PACK

6+ HOURS

BEFORE USE

Essential for optimal cooling



PRESS & HOLD 3s

Switch between
°C ↔ °F



1 Charge ≈ 1000 Uses

Up to 2 Months Battery Life*

*Battery life may vary based on usage.

⚠ Avoid direct sunlight, Keep in a cool place.



Made In China



WE' RE HERE FOR YOU

Scan to get the most from your product.



User Guide

Tips & best practices



Support

Get help fast



Warranty

Easy & hassle-free



EXTEND YOUR WARRANTY

GET +1 YEAR FREE!

Scan the QR code to register within
30 days of purchase.



Response within 24 hours



Diabetes_Care@outlook.com